

Product Datasheet

Recombinant Human High Mobility Group Protein B1/HMGB1 (N-MARI) (orb1648469)

| | |
|------------------------|---|
| Description | Recombinant Human High Mobility Group Protein B1/HMGB1 (N-MARI) |
| Species/Host | E. coli |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |