

Product Datasheet

INF2 antibody (orb164847)

Description	Rabbit polyclonal antibody to INF2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	INF2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	135624
Uniprot ID	Q27J81
NCBI	NP_116103.1 , NP_071934.3 , NP_001026884.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.