

Product Datasheet

Recombinant Human Ubiquitin-Conjugating Enzyme E2 D1/UBE2D1/UbcH5a (N-GST) (orb1648479)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 D1/UBE2D1/UbcH5a (N-GST)
Species/Host	E. coli
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.