

Product Datasheet

Recombinant Human Ubiquitin-Conjugating Enzyme E2 Variant 2/UBE2V2/DDVIT1 (N-6His) (orb1648582)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 Variant 2/UBE2V2/DDVIT1 (N-6His)
Species/Host	E. coli
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.