

Product Datasheet

Recombinant Human Ubiquitin-Conjugating Enzyme E2 C/UBE2C/UBCH10 (N-6His, T7 tag) (orb1648640)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 C/UBE2C/UBCH10 (N-6His, T7 tag)
Species/Host	E. coli
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.