

Product Datasheet

Recombinant Human Amphoterin-Induced Protein 2/AMIGO2/Alivin-1 (C-6His) (orb1649168)

| | |
|------------------------|---|
| Description | Recombinant Human Amphoterin-Induced Protein 2/AMIGO2/Alivin-1 (C-6His) |
| Species/Host | Human |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |