

## Product Datasheet

# Recombinant Human Platelet-Derived Growth Factor BB/PDGF-BB (orb1649481)

<b>Description</b>	Recombinant Human Platelet-Derived Growth Factor BB/PDGF-BB
<b>Species/Host</b>	E. coli
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.