

## Product Datasheet

# Recombinant Human Ubiquitin-Conjugating Enzyme E2 T/UBE2T (N-6His) (orb1649493)

<b>Description</b>	Recombinant Human Ubiquitin-Conjugating Enzyme E2 T/UBE2T (N-6His)
<b>Species/Host</b>	E. coli
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.