

Product Datasheet

Recombinant Human Fatty Acid-Binding Protein 4/FABP4/A-FABP/aP2 (N-6His) (orb1649540)

Description	Recombinant Human Fatty Acid-Binding Protein 4/FABP4/A-FABP/aP2 (N-6His)
Species/Host	E. coli
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Note	For research use only
Expiration Date	6 months from date of receipt.