

## Product Datasheet

# Recombinant Human Tumor Necrosis Factor- $\alpha$ /TNF- $\alpha$ High Active Mutant (orb1649618)

<b>Description</b>	Recombinant Human Tumor Necrosis Factor- $\alpha$ /TNF- $\alpha$ High Active Mutant
<b>Species/Host</b>	E. coli
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.