



Product Datasheet

CALCRL antibody (orb165003)


Description	Rabbit polyclonal antibody to CALCRL
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CALCRL
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52979
Uniprot ID	Q16602
NCBI	NP_005786.1 , NP_001258680.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot:Anti-CALCRL Antibody (Cente...



Western blot:Anti-CALCRL Antibody (Cente...



Immunohistochemical analysis of human he...