

Product Datasheet

Recombinant Human T-cell immunoreceptor with Ig and ITIM domains(TIGIT),partial (orb1650066)

| | |
|------------------------|---|
| Description | Recombinant Human T-cell immunoreceptor with Ig and ITIM domains(TIGIT),partial |
| Species/Host | E. coli |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | N-terminal 6xHis-SUMO-tagged |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |