

Product Datasheet

Recombinant Human Tumor necrosis factor receptor superfamily member 11B protein(TNFRSF11B) (Active) (orb1650511)

Description	Recombinant Human Tumor necrosis factor receptor superfamily member 11B protein(TNFRSF11B) (Active)
Species/Host	Yeast
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal FC-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.