

Product Datasheet

Recombinant Human Tumor necrosis factor receptor superfamily member 1B protein(TNFRSF1B) (Active) (orb1650516)

| | |
|------------------------|---|
| Description | Recombinant Human Tumor necrosis factor receptor superfamily member 1B protein(TNFRSF1B) (Active) |
| Species/Host | E. coli |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | NO-tagged |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |