

## Product Datasheet

# Recombinant Human Tumor necrosis factor receptor superfamily member 17 protein(TNFRSF17) (Active) (orb1650518)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human Tumor necrosis factor receptor superfamily member 17 protein(TNFRSF17) (Active)   |
| <b>Species/Host</b>    | E. coli   |
| <b>Storage</b>         | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Tag</b>             | NO-tagged   |
| <b>Note</b>            | For research use only   |
| <b>Expiration Date</b> | 6 months from date of receipt.  |