

## Product Datasheet

# Recombinant Rat Chemokine (C-X-C motif) ligand 12 (Stromal cell-derived factor 1) protein(Cxcl12) (Active) (orb1650559)

<b>Description</b>	Recombinant Rat Chemokine (C-X-C motif) ligand 12 (Stromal cell-derived factor 1) protein(Cxcl12) (A...
<b>Species/Host</b>	E. coli
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	NO-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.