

Product Datasheet

Recombinant Rat Chemokine (C-X-C motif) ligand 12 (Stromal cell-derived factor 1) protein(Cxcl12) (Active) (orb1650560)

Description	Recombinant Rat Chemokine (C-X-C motif) ligand 12 (Stromal cell-derived factor 1) protein(Cxcl12) (A...
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	NO-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.