

Product Datasheet

Recombinant Human Parathyroid hormone protein(PTH),15N Stable Isotope Labeled(Active) (orb1650572)

Description	Recombinant Human Parathyroid hormone protein(PTH),15N Stable Isotope Labeled(Active)
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	NO-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.