

Product Datasheet

Recombinant Rhesus Macaque Pro-interleukin-16 [Cleaved into: Interleukin-16 protein(IL16) (Active) (orb1650666)]

Description	Recombinant Rhesus Macaque Pro-interleukin-16 [Cleaved into: Interleukin-16 protein(IL16) (Active)]
Species/Host	E. coli
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	NO-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.