

## Product Datasheet

# Recombinant Human C-C motif chemokine 3 protein(CCL3) (Active) (orb1650799)

<b>Description</b>	Recombinant Human C-C motif chemokine 3 protein(CCL3) (Active)
<b>Species/Host</b>	E. coli
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	NO-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.