

Product Datasheet

Recombinant Human Myelin P2 protein(PMP2) (orb1651708)

| | |
|------------------------|---|
| Description | Recombinant Human Myelin P2 protein(PMP2) |
| Species/Host | E. coli |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Tag | N-terminal GST-tagged |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |