

## Product Datasheet

# Recombinant *Drosophila melanogaster* Ubiquitin-conjugating enzyme E2-17KDA(UbcD6) (orb1652488)

<b>Description</b>	Recombinant <i>Drosophila melanogaster</i> Ubiquitin-conjugating enzyme E2-17KDA(UbcD6)
<b>Species/Host</b>	<i>E. coli</i>
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Tag</b>	N-terminal 6xHis-SUMO-tagged
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.