

Product Datasheet

General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit (orb1653348)

Description	General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.118 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.