

## Product Datasheet

# Human BNIP3L / BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like ELISA Kit (orb1653435)

<b>Description</b>	Human BNIP3L / BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.14 ng/mL. Assay Type: Sandwich.
<b>Expiration Date</b>	Please enquire.