

Product Datasheet

Human Mucosal addressin cell adhesion molecule-1 (MADCAM-1) ELISA Kit (orb1653519)

Description	Human Mucosal addressin cell adhesion molecule-1 (MADCAM-1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 37 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.