

## Product Datasheet

# Human Ragulator complex protein LAMTOR3 (MAPKSP1) ELISA Kit (orb1653599)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Ragulator complex protein LAMTOR3 (MAPKSP1) ELISA Kit  |
| <b>Reactivity</b>          | Human  |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |