

Product Datasheet

Goat Malondialdehyde (MDA) ELISA Kit (orb1653704)

Description Goat Malondialdehyde (MDA) ELISA Kit

Reactivity Goat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: Request Information.
Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.