

Product Datasheet

Human DNA-3-methyladenine glycosylase (MPG) ELISA Kit (orb1654057)

Description	Human DNA-3-methyladenine glycosylase (MPG) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.