

Product Datasheet

Sheep NADH-ubiquinone oxidoreductase chain 1 (MT-ND1) ELISA Kit (orb1654295)

Description	Sheep NADH-ubiquinone oxidoreductase chain 1 (MT-ND1) ELISA Kit
Reactivity	Sheep
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.