

Product Datasheet

Human Melatonin receptor type 1A (MTNR1A) ELISA Kit (orb1654355)

Description	Human Melatonin receptor type 1A (MTNR1A) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.