

Product Datasheet

Human Prenylcysteine oxidase 1 (PCYOX1) ELISA Kit (orb1654865)

Description	Human Prenylcysteine oxidase 1 (PCYOX1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.118 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.