

Product Datasheet

RRBP1 antibody (orb165492)

Description	Rabbit polyclonal antibody to RRBP1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RRBP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	152456
Uniprot ID	Q9P2E9
NCBI	NP_001036041.1 , NP_004578.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

