

## Product Datasheet

# Rat Peroxiredoxin-6 (PRDX6) ELISA Kit (orb1655166)

**Description** Rat Peroxiredoxin-6 (PRDX6) ELISA Kit

**Reactivity** Rat

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: Request Information.  
Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200  $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.

**Expiration Date** Please enquire.