

Product Datasheet

Human Soluble endothelial protein C receptor (sEPCR) ELISA Kit (orb1655315)

Description	Human Soluble endothelial protein C receptor (sEPCR) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 ng/mL. Sensitivity: 1.95 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.