

Product Datasheet

Human Protein prune homolog (PRUNE) ELISA Kit (orb1655415)

Description	Human Protein prune homolog (PRUNE) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 47-3000 pg/mL. Sensitivity: 11.75 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.