

Product Datasheet

Mouse Parathyroid Hormone Related Protein (PTHrP) ELISA Kit (orb1655494)

Description	Mouse Parathyroid Hormone Related Protein (PTHrP) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.2 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.