

Product Datasheet

Human Protein tyrosine phosphatase 1B (PTP1B) ELISA Kit (orb1655518)

Description	Human Protein tyrosine phosphatase 1B (PTP1B) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.054 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.