

## Product Datasheet

### Human Putative peptide YY-3 (PYY3) ELISA Kit (orb1655615)

|                            |                                                                                                                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Putative peptide YY-3 (PYY3)<br>ELISA Kit                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                      |
| <b>Note</b>                | For research use only                                                                                                                                                                                      |
| <b>Application notes</b>   | Detect Range: 0.156-10 pg/mL.<br>Sensitivity: 0.078 pg/mL. Assay<br>Time: 1-4.5 h. Detection<br>wavelength: 450 nm. Sample<br>Volume: 1-200 µL. Sample Type:<br>Serum, Plasma, Other biological<br>fluids. |
| <b>Expiration Date</b>     | Please enquire.                                                                                                                                                                                            |