

Product Datasheet

Rat Acetylcholine (ACH) ELISA Kit (orb1655831)

| | |
|----------------------------|--|
| Description | Rat Acetylcholine (ACH) ELISA Kit |
| Reactivity | Rat |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 12.35-1000 pg/mL. Sensitivity: 5.29 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |