

Product Datasheet

Human Canalicular multispecific organic anion transporter 1 (ABCC2) ELISA Kit (orb1655864)

| | |
|----------------------------|--|
| Description | Human Canalicular multispecific organic anion transporter 1 (ABCC2) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.156-10 ng/mL. Sensitivity: 0.067 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |