

## Product Datasheet

# Mouse Alpha-2-macroglobulin (A2M) ELISA Kit (orb1655886)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse Alpha-2-macroglobulin (A2M)<br>ELISA Kit  |
| <b>Reactivity</b>          | Mouse   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 3.70-300 ng/mL.<br>Sensitivity: 1.51 ng/mL. Assay Time:<br>1-4.5 h. Detection wavelength: 450<br>nm. Sample Volume: 1-200 $\mu$ L.<br>Sample Type: Serum, Plasma, Other<br>biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |