

Product Datasheet

Rat 8-Hydroxy-desoxyguanosine (8-OHdG) ELISA Kit (orb1655889)

Description	Rat 8-Hydroxy-desoxyguanosine (8-OHdG) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.078 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.