

Product Datasheet

Human 8-iso prostaglandin F2 α (8-iso-PGF2a) ELISA Kit (orb1655893)

Description	Human 8-iso prostaglandin F2 α (8-iso-PGF2a) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.25-2000 pg/mL. Sensitivity: 15 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.