

Product Datasheet

Human 25-hydroxy vitamin D3 (25 HVD3) ELISA Kit (orb1655910)

Description	Human 25-hydroxy vitamin D3 (25 HVD3) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.12-200 ng/mL. Sensitivity: 1.56 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.