

## Product Datasheet

### Rat Relaxin 1 (RLN1) ELISA Kit (orb1655968)

<b>Description</b>	Rat Relaxin 1 (RLN1) ELISA Kit
<b>Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 2.4-150 pg/mL. Sensitivity: 0.6 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.
<b>Expiration Date</b>	Please enquire.