

Product Datasheet

Human Amiloride-sensitive sodium channel subunit alpha (SCNN1A) ELISA Kit (orb1656043)

Description	Human Amiloride-sensitive sodium channel subunit alpha (SCNN1A) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.