

Product Datasheet

Mouse Semaphorin-4D (SEMA4D) ELISA Kit (orb1656083)

| | |
|----------------------------|--|
| Description | Mouse Semaphorin-4D (SEMA4D) ELISA Kit |
| Reactivity | Mouse |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.78-50 ng/mL. Sensitivity: 0.30 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |