

Product Datasheet

Human Sialic acid-binding Ig-like lectin 6 (SIGLEC6) ELISA Kit (orb1656143)

Description	Human Sialic acid-binding Ig-like lectin 6 (SIGLEC6) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.057 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.