

Product Datasheet

Rat Serotonin transporter (5-HTT/SERT) ELISA Kit (orb1656219)

Description	Rat Serotonin transporter (5-HTT/SERT) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.061 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.