

## Product Datasheet

# Human Soluble lectin-like oxidized low density lipoprotein receptor-1 (sLOX-1) ELISA Kit (orb1656244)

<b>Description</b>	Human Soluble lectin-like oxidized low density lipoprotein receptor-1 (sLOX-1) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 62.5-4000 pg/mL. Sensitivity: 38 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
<b>Expiration Date</b>	Please enquire.